

How Do I Get Clear Autoclave Bags?

Buy clear autoclave bags for your lab by ordering online

or

Pick them up in the Fisher Stockroom (LGRT):

Item Name	Fisher #	#/Box	Size
Autoclave bags/clear	01-814-2	200	24" x 30"
Autoclave bags/clear	01-814-3	200	24" x 36"



Use only *autoclave* bags in an autoclave. Other bags melt.

Do not use red or orange autoclave bags.

Biohazard labels are not allowed.

Questions? Contact: askehs@umass.edu